REFERENCES

- Polson RJ; Kenna JG; Shears IP; Bomford A; Williams R. Measurement of ferritin in serum by an indirect competitive enzyme-linked immunosorbant assay. Clin Chem 1988; 34(4):661-4.
- Iacobello C; Ruggeri G; Signorini C; Di Cello GC; Albertini A. A semi-automated ELISA for serum ferritin determination. Clin Chem 1987; 32(5):899.
- Aulbert E; Steffens O. Serum ferritin--a tumor marker in malignant lymphomas? Onkologie 1990; 13(2):102-8.
- Essen A; Ozen H; Ayhan A; Ergen A; Tasar C; Remzi F. Serum ferritin: a tumor marker for renal cell carcinoma. J Urol 1990;145(6):1134-7.
- Borque L; Rus A; del Cura J; Maside C; Escanero J. Automated quantitative nephelometric latex immunoassay for determining ferritin in human serum. J Clin Lab Anal 1992; 6(4):239-44
- 6. Dawson DW; Fish DI; Shackleton P The accuracy and clinical interpretation of serum ferritin assays. Clin Lab Haematol 1992; 14(1):47-52.

2025-08-08



INTENDED USE

The Ferritin ELISA Kit is intended for the quantitative measurement of Ferritin in human serum or plasma. For Research Use Only. For professional use only. Not for use in diagnostic procedures.

SUMMARY AND EXPLANATION

Ferritin is a large molecule with a molecular weight of approximately 450,000 Daltons, and consists of a protein shell around an iron core; each molecule of Ferritin may contain as many as 4,000 iron atoms. Ferritin is the major intracellular iron-storage protein. Research applications include iron homeostasis, oxidative stress, inflammation-mediated iron redistribution, and ferritin gene regulation studies.

PRINCIPLE OF THE TEST

This Ferritin ELISA kit is a solid phase sandwich assay method, based on a streptavidin-biotin principle. The standards, samples and the biotinylated Anti-Ferritin antibody reagent are added into designated wells, coated with Streptavidin. Endogenous Ferritin in the patient's serum binds to the antigenic site of the biotinylated Anti-Ferritin antibody. Simultaneously, the biotinylated antibody is immobilized onto the wells through the high affinity Streptavidin-Biotin interaction. Unbound protein and excess biotin conjugated antibody are washed off by wash buffer. Upon the addition of the Peroxidase (HRP) conjugated Anti-Ferritin antibody reagent, a sandwich complex is formed, the analyte of interest being in between the two highly specific antibodies, labeled with Biotin and HRP. Unbound protein excess enzyme conjugated antibody reagent is washed off by wash buffer. Upon the addition of the substrate, the intensity of color developed is directly proportional to the concentration of Ferritin in the samples. A standard curve is prepared relating color intensity to the concentration of the Ferritin.

	MATERIALS PROVIDED	96 TESTS
1.	Microwells coated with Streptavidin	12x8x1
2.	Ferritin Standards: 6 vials (ready to use)	0.5 mL
3.	Ferritin Controls: 2 vials (ready to use)	0.5 mL
4.	Anti-Ferritin Biotin Reagent (ready to use)	12 mL
5.	Anti-Ferritin Enzyme Reagent (ready to use)	12 mL
6.	TMB Substrate: 1 bottle (ready to use)	12 mL
7.	Stop Solution: 1 bottle (ready to use)	12 mL
8.	20X Wash Concentrate: 1 bottle	25 mL

MATERIALS NOT PROVIDED

- 1. Distilled or deionized water
- 2. Precision pipettes
- 3. Disposable pipette tips
- 4. ELISA reader capable of reading absorbance at 450nm
- Absorbance paper or paper towel
- 6. Graph paper

STORAGE AND STABILITY

- 1. Store the kit at 2 8° C.
- 2. Keep microwells sealed in a dry bag with desiccants.
- 3. The reagents are stable until expiration of the kit.
- 4. Do not expose reagent to heat, sun, or strong light.



Calbiotech, Inc.

1935 Cordell Ct., El Cajon, CA 92020 USA Tel (619) 660-6162 | Fax (619) 660-6970 | Web www.calbiotech.com

WARNINGS AND PRECAUTIONS

Potential biohazardous materials:

- The standards contain human source components which have been tested and found non-reactive for hepatitis B surface antigen as well as HIV antibody with FDA licensed reagents.
 However, as there is no test method that can offer complete assurance that HIV, Hepatitis B virus or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2, as recommended in the Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories." 1984.
- Do not pipette by mouth. Do not smoke, eat, or drink in the areas in which specimens or kit reagents are handled.
- The components in this kit are intended for use as an integral unit. The components of different lots should not be mixed.
- 4. It is recommended that standards, control and serum samples be run in duplicate
- Optimal results will be obtained by strict adherence to this protocol. Accurate and precise
 pipetting, as well as following the exact time and temperature requirements prescribed are
 essential. Any deviation from this may yield invalid data.

SPECIMEN COLLECTION HANDLING

This assay is designed for use with human serum or plasma samples obtained in compliance with applicable laws, regulations, and institutional policies. Handle and store samples using procedures appropriate for research use. Samples may be stored refrigerated (2-8 °C) for up to seven days, or frozen (-20 °C or below) for up to six months. Avoid repetitive freeze—thaw cycles.

REAGENTS PREPARATION

Prepare 1X Wash buffer by adding the contents of the bottle (25 ml, 20X) to 475 ml of distilled or deionized water. Store at room temperature (20-25°C).

ASSAY PROCEDURE

Prior to assay, allow reagents to stand at room temperature (20-25°C). Gently mix all reagents before use.

- 1. Place the desired number of coated strips into the holder
- 2. Pipette 25 µl of Ferritin standards, controls, and samples in to appropriate wells.
- 3. Add 100 μ l of Biotin Reagent into each well. Shake the plate for (10-30) sec.
- 4. Cover the plate and incubate for 30 minutes at room temperature (20-25° C).
- 5. Remove liquid from all wells. Remove liquid from all wells. Wash wells three times with 300 μ l of 1X wash buffer. Blot on absorbance paper or paper towel.
- 6. Add 100 μl of Enzyme Reagent into each well.
- 7. Cover the plate and incubate for 30 minutes at room temperature (20-25° C).
- 8. Remove liquid from all wells. Remove liquid from all wells. Wash wells three times with 300 μ l of 1X wash buffer. Blot on absorbance paper or paper towel.
- Add 100 μl of TMB substrate to all wells.
- 10. Incubate for 15 minutes at room temperature.
- 11. Add 50 µl of stop solution to all wells. Shake the plate 10-20 seconds to mix the solution.
- Read absorbance on ELISA Reader at 450 nm within 15 minutes after adding the stopping solution.

CALCULATION OF RESULTS

The standard curve is constructed as follows:

- Check Ferritin standard value on each standard vial. This value may vary from lot to lot. Make sure you check the value on every kit. See example of the standard below.
- To construct the standard curve, plot the absorbance for the Ferritin standards (vertical axis) versus the Ferritin standard concentrations (horizontal axis) on a linear graph paper. Draw the best curve through the points.
- 3. Read the absorbance for controls and each unknown sample from the curve. Record the value for each control or unknown sample.

IFU-FR248T-RC-V2

Example of a Standard Data

	OD 450 nm	Conc. ng/mL
Std 1	0.019	0
Std 2	0.144	10
Std 3	0.531	50
Std 4	1.197	150
Std 5	1.987	400
Std 6	2.543	800

LIMITATIONS OF THE TEST

1. Do not use sodium azide as preservative. Sodium azide inhibits HRP enzyme activities.